

Traceable® Water-/Shock-Resistant Stopwatch model 1045



Specifications: Timing range is 23 hours, 59 minutes, 59 seconds · Resolution is 1-second · Accuracy is 0.01% · Size is 2-1/8 x 2-3/8 x 1/2 inches · Weight is 1 ounce

General-purpose LCD Digital Stopwatch times to 24 hours. Features include: single-action, time-out, cumulative-split, and continuous timing. Timing resolution is 1/100 of a second for the first 30 minutes, then continues timing to 23 hours, 59 minutes, 59 seconds, with a 1-second resolution.

As a chronograph, the unit displays time of day in hours, minutes, and seconds, A.M./P.M., date, and day of the week. Also has a time-of-day alarm. Tactilefeel, "click" switches provide a positive action. Accuracy is 0.01%. Digits are 1/4-inch high. Stopwatch is 2-1/8 x 2-3/8 x 1/2 inches and weight is 1 ounce.

To assure accuracy an individually serial-numbered Traceable® Certificate is provided from our ISO 17025 calibration laboratory accredited by A2LA. It indicates traceability to standards provided by NIST (National Institute of Standards and Technology). Supplied also with, 2-year "always-on" battery and lanyard.