



Visual-Audible Bench Timer

- Entire LCD illuminates and flashes on and off, hard to miss even in a noisy environment.
- Loud audible alarm.
- Large digits & off button.
- Sturdy rubberized case.



Scan this QR code to view our product demo.



No.	Description	Each
810014	Visual-Audible Bench Timer	
810014C	Visual-Audible Bench Timer Certified	



Digital Hourglass

- The simplicity of an hourglass adjustable to any interval.
- Shows min and sec remaining at the top of the display, and min and sec elapsed at the bottom.
- Restarted simply by flipping the timer over.



No.	Description	Each
810044	Digital Hourglass	
810044C	Digital Hourglass Certified	



Mini Timer

- Compact and easy-to-hold.
- Touch Tone buttons.
- Magnetic back.



No.	Description	Each
810019	Mini Timer	
810019C	Mini Timer Certified	



Stopwatch Timer

- True stopwatch mode with 1/100 second resolution.
- Lap and split.
- Timer mode.



No.	Description	Each
810025R	Stopwatch Timer	
810025RC	Stopwatch Timer Certified	



White Large Display Timer

- Easily read from across a room
- Simultaneously displays either one count-up channel and a clock, or two count-down channels.
- Extra large start/stop button and a loud alarm.
- Wall mount.



No.	Description	Each
810005	White Large Display Timer	
810005C	White Large Display Timer Certified	



Project Timer

- Counts down to completion of a project up to 8 years in the future in days, hours, minutes, seconds, and 1/100 seconds.
- View either the target date or the actual count down.
- Insert the name of your project into a special window.



No.	Description	Each
810017	Project Timer	
810017C	Project Timer Certified	



Count Up/Count Down

- Numeric keypad allows for fast programming.
- Counts up or down for 99 minutes (810024R) or 100 hours (810026R) in intervals of 1 second.



No.	Description	Each
810024R	Count Up/Count Down 99 Min	
810024RC	Count Up/Count Down 99 Min Certified	
810026R	Count Up/Count Down 100 Hr	
810026RC	Count Up/Count Down 100 Hr Certified	



Simultaneous Start Timer

- 4 channels can be set independently but started simultaneously.
- Loud alarm.



No.	Description	Each
810027	Simultaneous Start Timer	
810027C	Simultaneous Start Timer Certified	



Visual Alarm Timer

- Select a bright flashing yellow LED or an audible alarm.
- Counts time elapsed after the alarm.



No.	Description	Each
810040	Visual Alarm Timer 99 Min	
810040C	Visual Alarm Timer 99 Min Certified	
810041	Visual Alarm Timer 99 Hr	
810041C	Visual Alarm Timer 99 Hr Certified	



5 Channel Timer

- Large Display.
- Indicates time elapsed after countdown reaches zero.
- Has a stand, pocket clip and magnet.



No.	Description	Each
810015	5 Channel Timer	
810015C	5 Channel Timer Certified	



Pocket Large Display Timer

- ¾ inch high digits.
- Easily portable.
- Warning alarms sound at 3 and 5 minutes prior.
- Records elapsed time following the alarm.



No.	Description	Each
810046	Pocket Large Display Timer	
810046C	Pocket Large Display Timer Certified	



Continuous Alarm Timer

- A bright, flashing yellow LED visual alarm, or an audible alarm.
- Will flash or sound continuously for weeks until shut off.
- Counts time elapsed after the alarm.



No.	Description	Each
810045	Continuous Alarm Timer	
810045C	Continuous Alarm Timer Certified	



Talking Timer

- Eliminate the need to continually visually check time remaining.
- Audible count-down can be heard throughout the lab.
- Useful for work in low light.



No.	Description	Each
810038	Talking Timer	
810038C	Talking Timer Certified	

Model:	Count Up Channel	Count Down Channel	Range	Resolution	Recalls Previous Setting	Clock	Alarm Clock
Visual-Audible Bench Timer 810014	1	1	99 Min / 59 Sec	1 Sec	✓		
White Large Display Timer 810005	1	2	24 Hrs	1 Sec	✓	✓	
Digital Hourglass 810044	0	1	100 Min	1 Sec	✓		
Project Timer 810017	1	2	24Hrs/Clock up to 8 years	1/100 Sec		✓	
Mini Timer 810019	0	1	99 Min	1 Sec	✓		
Count Up/Count Down 810024R	1	1	99 Min	1 Sec	✓		
Count Up/Count Down 810026R	1	1	100 Hrs	1 Sec	✓	✓	✓
Stopwatch Timer 810025R	1	1	100 Hrs	1/100 Sec	✓		
Simultaneous Start Timer 810027	4	4	100 Hrs	1 Sec	✓	✓	✓
	1	0	10 Hrs	1/100 Sec			
Pocket Large Display Timer 810046	1	1	100 Min	1 Sec	✓		
Visual Alarm Timer 810040:	0	1	99 Min 59 Sec	1 Sec	✓		
810041:	0	1	99 Hours 59 Min	1 Sec			
Continuous Alarm Timer 810045	1	1	9999 Min	1 Min / 1 Sec	✓		
5 Channel Timer 810015	5	4	100 Hrs	1 Sec	✓	✓	
Talking Timer 810038	1	1	24 Hrs	1 Sec	✓	✓	