

3000 Series Dry-Well



Portable Mini Heat Source

Model #271-3XX

0.1° resolution

Small, lightweight heat source

Quick tester to 572°F


Accepts all standard probe diameters from 1/8" to 3/8" dia.

Four block configurations to choose from


NIST-Traceable standard

thermometer


3001 Block Configuration




3002 Block Configuration



3003 Block Configuration




3004 Block Configuration



Removable Brass Inserts for 3004

271-321



271-322



271-323



271-324



271-325



Specification

Range	91.4°F to 572°F (33°C to 300°C)
Accuracy	±0.9°F (91.4 to 391.8°F) & ± 1.8°F (392 to 572°F) ±0.5°C (33 to 199°C) & ± 1°C (200 to 300°C)
Resolution	0.1°
Stability	±0.5°C
Heating Time	Ambient to 572°F (300°C) after 10 minutes
Cooling Time	572°F (300°C) to 212°F (100°C) after 10 minutes
Stabilization	5 minutes
Well Depth	4" (100mm)
Power	230/115 volts
Display	0.4" (10mm) 4-digit LED
Dimensions	2.25" x 5" x 6.25" (57 x 125 x 158mm)
Weight	2 lbs. (900 g.)

